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
Clinical Update

Number Thirteen

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1. Fischbach FT. A Manual of laboratory and Diagnostic Tests, 7th Ed. Lippincott, Williams & Wilkins, Philadelphia 2004:404-05
2. Corbet JV. Laboratory Tests and Diagnostic Procedures, 6th Ed. Pearson Prentice Hall, Upper Saddle River, New Jersey 2004:286-87
3. Morishita Y, Arikawa K, Yamashita M, et al. Fatal hemolysis due to unidentified causes following mitral valve replacement with bileaflet tilting disc valve prosthesis. Heart Vessels 1987;3:100-03
4. Barmada H, Starr A. Clinical hemolysis with the St. Jude valve without paravalvular leak. Medical Progress Through Technology 1994;20:191-94
5. Skoulargis J, Essop MR, Skudicky D, et al. Frequency and severity of intravascular hemolysis after left-sided cardiac valve replacement with Medtronic Hall and St. Jude Medical prostheses, and influence of prosthetic type, position, size and number. Am J Cardiol 1993;71:587-91
6. Mecozzi G, Milano AD, De Carlo M, et al. Intravascular hemolysis in patients with new-generation prosthetic heart valves: a prospective study. J Thorac Cardiovasc Surg 2002;123:550-56
7. Lund O, Emmersten K, Nielsen TT, et al. Impact of size mismatch and left ventricular function on performance of the St. Jude disc valve after aortic valve replacement. Ann Thorac Surg 1997;63:1227-34
8. Niinami H, Aomi S, Tomioka H, et al. A comparison of the in vivo performance of the 19-mm St. Jude medical hemodynamic plus and 21-mm standard valve. Ann Thorac Surg 2002;74:1120-24
9. Ismeno G, Renzulli A, De Feo M, et al. Standard versus hemodynamic plus 19-mm St. Jude Medical aortic valves. J Thorac Cardiovasc Surg 2001;121:723-28
10. Birnbaum D, Laczkovics A, Heidt M, et al. Examination of hemolytic potential with the On-X® prosthetic heart valve. J Heart Valve Dis 2000;9:142-45
11. Summary of Safety and Effectiveness, On-X® prosthetic heart valve, FDA PMA P000037 and P000037/S1, May 30, 2001 and March 6, 2002, European Primary Trial, updated May 31, 2003
12. Ely J, Emken M, Accuntius J, et al. Pure pyrolytic carbon; preparation and properties of a new material, On-X carbon, for mechanical heart valve prostheses. J Heart Valve Dis 1998; 7: 626-32
13. Bokros JC, Ely JL, Emken MR, et al. Prosthetic heart valve with improved blood flow. US Patent No. 5,772,694: Issued June 30, 1998
14. Bokros JC, Emken MR, Haubold AD et al. Prosthetic heart valve. US Patent No. 5,308,361: Issued May 3, 1994 and US Patent No. 5,137,532: Issued Aug. 11, 1992
15. Bokros JC, Stupka JC, More RB. Prosthetic heart valve. US Patent No. 6,096,075: Issued Aug. 1, 2003
16. Bokros JC, Ely JL, Emken MR et al. Prosthetic heart valve. US Patent No. 5,642,324: Issued June 24, 1997
17. Gencbay M, Degertekin M, Basaran Y, et al. Microbubbles associated with mechanical heart valves: their relation with serum lactic dehydrogenase levels. Am Heart J 1999;137:463-68
18. Rambod E, Gharib M, Denis JL, et al. Can the closure impact of mechanical heart valves (MHV) release microbubbles in gas-saturated solutions? An in vitro study. J Am Coll Cardiol 1996;27(suppl A):234A
19. Williams MA, Crause L, van Riet S. A comparison of mechanical valve performance in a poorly anticoagulated community. J Card Surg 2004;19 (5):410-14
20. Moidl R, Simon P, Wolner E, et al. The On-X® prosthetic heart valve at five years. Ann Thorac Surg 2002;74:S1312-17
21. Laczkovics A, Heidt M, et al. Early clinical experience with the On-X® prosthetic heart valve. J Heart Valve Dis 2001;10:94-99
22. Cavitation Resistance Chart. Austin, Texas 2000: Medical Carbon Research Institute. Data generated by NHC Wang, University of Miami and H Reul, Helmholtz Institute, Aachen, Germany.

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ON-X Life Technologies, Inc.