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Clinical Update

Number Fourteen

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1. Bonow RO, Carabello B, de leon AC, et al. ACC/AHA Pocket Guideline for Management of Patients with Valvular Heart Disease. July 2000. Adapted from JACC 1998; 32(5):1486-1588
2. Ghosh P. Changing trends in valve surgery in Europe: 1991-2000. J Heart Valve Dis 2004; 13:1-10
3. Emery RW, Krogh CC, Arom KV, et al. The St. Jude Medical cardiac valve prosthesis: a 25-year experience with single valve replacement. Ann Thorac Surg 2005; 79: 776-83
4. Rahimtoola SH. Choice of prosthetic heart valve for adult patients. JACC 2003;41 (6):893-904.
5. Grunkemeier GL, Li H-H, Naftel DC, et al. Long-term performance of heart valve prosthesis. Curr Probl Cardiol 2000;25:73-156
6. Savit M, Gao G, Wu Y, Grunkemeier GL, et al. Twelve-year comparison between the Carpentier-Edwards pericardial and bileaflet mechanical valves in the aortic position. Poster presentation at the 41st Annual Meeting of the Society of Thoracic Surgeons 2005;P-344
7. Hammermeister KE, Sethi GK, Henderson WC, et al. Outcomes 15 years after valve replacement with a mechanical versus bioprosthetic valve: final report of the Veterans Affairs randomized trial. JACC 2000;36:1152-58
8. Bloomfield P, Wheatley DJ, Prescott RJ, et al. Twelve-year comparison of a Bjork-Shiley mechanical heart-valve with porcine bioprosthesis. N Engl J Med 1991; 324: 573-79
9. Philips SJ. Searching for the truth: a mechanical or a tissue valve? J Heart Valve Dis 2004;13(Suppl.1):S95-S98
10. Wheatley D. The "Threshold Age" in choosing biological versus mechanical prostheses in western countries. J Heart Valve Dis 2004;13 (Suppl.1):S91-S94
11. Piehler J, Blackstone EH, Bailey KR, et al. Reoperation on prosthetic heart valve-patient specific estimation of in-hospital events. J Thorac Cardiovasc Surg 1995;109:30-48
12. Kirsch M, Nakashima K, Kubota S, et al. The risk of reoperative heart valve procedures in octogenarian patients. J Heart Valve Dis 2004;13:991-96
13. Moidl R, Simon P, Wolner E, et al. The On-X® Prosthetic Heart Valve at five years. Ann Thorac Surg 2002;74:S1312-17
14. Williams MA, Crause L, van Riet S. A comparison of mechanical valve performance in a poorly anticoagulated community. J Card Surg 2004;19 (5):410-14
15. Laczkovics A, Heidt M, Oelert H, et al. Early clinical experience with the On-X® Prosthetic Heart Valve. J Heart Valve Dis 2001;10:94-99
16. Kortke H, Korfer R. International normalized ratio self-management after mechanical heart valve replacement: is an early start advantageous? Ann Thorac Surg 2001;72:44-48
17. Summary of Safety and Effectiveness, On-X® Prosthetic Heart Valve, FDA PMA P000037, May 30,2000, updated through May 31, 2003
18. Dellgren G, David TE, Raanani E, et al. Late hemodynamic and clinical outcomes of aortic valve replacement with the Carpentier-Edwards Perimount pericardial bioprosthesis. J Thorac Cardiovasc Surg 2002;124:146-54
19. Marchand MA, Aupart MR, Norton R, et al. Fifteen-year experience with the mitral Carpentier-Edwards PERIMOUNT pericardial bioprosthesis. Ann Thorac Surg 2001; 71:S236-39
20. Cartier PC, Dumesnil JG, Metras J, et al. Clinical and hemodynamic performance of the Freestyle aortic root bioprostheses. Ann Thorac Surg 1999;67:345-51
21. Novoa R, Sanofsky S, Sandhu A. Prevalence of patient-prosthesis mismatch in a community hospital. Medical Carbon Research Institute, LLC, Austin, TX USA: Innovations in Cardiac Surgery 2003:27

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