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Clinical Update

Number Sixteen

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REFERENCES:

1. Wu ZJ, Gao BZ, Slonin JH, et al. Bileaflet mechanical heart valve at low cardiac output. ASAIOJ 1996;42(5):747-49
2. Geometric orifice areas(GOA) for all valve models with valve leaflets in place were calculated at Medical Carbon Research Institute, LLC, in Austin, Texas USA.
3. Shipkowitz T, Ambrus J, Kurk J, et al. Evaluation technique for bileaflet mechanical valves. J Heart Valve Dis 2002;11:275-82
4. Hwang NHC, Reul H, Reinhard P. In vitro evaluation of the long body On-X® bileaflet heart valve. J Heart Valve Dis 1998;7(5):561-68
5. Otto CM. Echocardiographic evaluation of prosthetic valve dysfunction. In Otto CM (ed), Textbook of Clinical Echocardiography, 2nd edition WB Saunders Co, Philadelphia, 2003:310
6. Klepetko W, Moritz A, Mlczech J, et al. Leaflet fracture in Edwards-Duromedics bileaflet valves. J Thorac Cardiovasc Surg 1989;97:90-94
7. Birnbaum D, Laczkovics A, Heidt M, et al. Examination of hemolytic potential with the On-X® Prosthetic Heart Valve. J Heart Valve Dis 2000;9:142-45
8. Gencbay M, Degertekin M, Basaran Y, et al. Microbubbles associated with mechanical heart valves: their relation with serum lactic dehydrogenase levels. Am Heart J 1999;137:463-68
9. Summary of Safety and Effectiveness, On-X® Prosthetic Heart Valve, FDA PMA P000037, May 30, 2000, updated through 2004
10. Kortke H, Korfer R. International normalized ratio self-management after mechanical heart valve replacement: is an early start advantageous? Ann Thorac Surg 2001;72:44-48
11. Summary of safety and effectiveness for the Carbomedics® Prosthetic Heart Valve, FDA PMA P900060, September 29, 1993
12. Chambers J, Laczkovics A, Greve H, et al. Early postoperative echocardiographic hemodynamic performance of the On-X® Prosthetic Heart Valve — a multi-center study. J Heart Valve Dis 1998;7(5):569-73
13. Björnerheim R, Ihlen H, Simonsen S, et al. Hemodynamic characterization of the Carbomedics mitral valve prosthesis. J Heart Valve Dis 1997;6:115-22
14. Aoyagi S, Yasunaga H, Sato T, et al. Doppler echocardiographic evaluation of the St. Jude Medical Valve. Artif Organs Today 1995;5(1):49-56